Biatain® Ag

Biatain® Plata

Non-Adhesive Foam Dressing with silver
Adhesive Foam Dressing with silver
Adhesive Foam Dressing with silver:
Sacral
Heel
Biatain Ag Foam Filler/Cavity with silver

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Manufacturer : Coloplast A/S 3050 DK
23311747 V3 The Coloplast logo is a registered trademark of Coloplast A/S.
© [2011-10.] All rights reserved Coloplast A/S, 3050 Humlebæk, Denmark.

Distributed by Coloplast Corp
Minneapolis, MN 55411, USA
Country of origin : Hungary

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Warning

Reuse of this single use product may create a potential risk to the user. Reprocessing, cleaning, disinfection and sterilization may compromise product characteristics which in turn create an additional risk of physical harm to or infection of the patient.

Precautions

A healthcare professional should frequently inspect and manage diabetic wounds and wounds, which are solely or partially caused by arterial insufficiency, in accordance with local standards. Biatain Ag foam dressings should not be used on patients with a known sensitivity to silver. In case of suspected allergic reaction, contact Coloplast A/S for further information. Do not use Biatain Ag foam dressings with oxidizing solutions e.g. hypochlorite and hydrogen peroxide solutions. Ensure that any other evaporating solution is completely dried off before dressing application. Biatain Ag foam dressings may cause a transient discolouration of the wound bed, which can be removed by gentle washing. Biatain Ag foam dressings must be removed prior to radiation treatment or examinations that include x-rays, ultrasonic treatment, diathermy or microwaves. Biatain Ag foam dressings are MR Safe and MR compatible up to 3 Tesla and can therefore be left in place during an MR scan (Tesla: The unit for the strength of a magnetic field). The use of Biatain Ag foam dressings during pregnancy, lactation and on children has not been demonstrated. Application of Biatain Ag foam dressings with enzymatic debriding agents has not been demonstrated. The use of Purlion Gel with Biatain Ag Cavity dressings has not been demonstrated. Biatain Ag foam cavity dressing is not for use in dry wounds or deep cavity wounds with only a small opening.

Do not use if package is damaged.

Do not reuse.

Biatain Ag Adhesive Foam Dressings should be stored in a horizontal position.

Biatain Ag foam dressings should be stored at room temperature.

Keep away from sunlight.

Coloplast accepts no liability for any injury or other loss that may arise in the use of this product in a manner contrary to Coloplast's current recommendations.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
Application

-  1. Biatain Ag Non-Adhesive Foam Dressings
The dressing must be applied with the plain (non-printed) side towards the wound.
-  2. Biatain Ag Adhesive Foam Dressings
Use the handles to prevent touching the dressing and to ensure aseptic application. Remove the protective film first and place the foam part of the dressing over the wound. Then remove the release films on the sides as part of the dressing application.
-  3. Biatain Ag Adhesive Foam Dressings / Sacral
Use the handles to ensure aseptic application. Remove the center protective film. Place the narrow end of the dressing as far down over the wound as possible and fix the dressing working up and outwards. Remove the handles.
-  4. Biatain Ag Adhesive Foam Dressings / Heel
Use the handles to ensure aseptic application. Remove the center protective film. The dressing is shaped like an "arrow". Fold the dressing in a 90° angle between the "arrow head" and "arrow tail". Let the "arrow" point away from the heel tip and fix the "arrow tail" first. Secondly, fix the "arrow head". Remove the handles one by one and fix the sides carefully, in order that the foam parts

touch or overlap.

5. Biatain Ag Foam Filler/Cavity Remove the paper backing. Loosely pack Biatain Ag cavity to fill approximately half the volume of the wound. On absorption of exudate Biatain Ag cavity expands to completely fill the cavity in a short period of time. Mold Biatain Ag cavity aseptically according to the wound shape by loosely folding the product along the perforations or by rolling it up.

To remove the dressing, gently lift the corners of the dressing away from the wound.

To remove the adhesive dressings, it is recommended to loosen the adhesive border before lifting the dressing away from the wound.

Biatain Ag sacral foam dressings should be removed from the top edge and down towards the anus to minimize the risk of transmitting infection.

Description

Biatain Ag foam dressings are sterile single-use, soft, highly absorbent and conformable antibacterial polyurethane foam dressings with ionic silver as the active component. Biatain Ag foam dressings provides an optimal moist wound healing environment and effective exudate management.

Biatain Ag	Dressing Size	Silver Content in mg.
Adhesive	5"x5"	60
Adhesive	7"x7"	160
Non-adh	4"x4"	95
Non-adh	6"x6"	215
Sacral	9"x9"	120
Heel	7"x8"	85

Biatain Ag foam dressings contain an antibacterial silver complex homogeneously dispersed throughout the foam. In the presence of exudate the release of ionic silver will continue for the recommended wear time up to 7 days.

Biatain Ag foam dressings are protected with semi-permeable film backings that are waterproof and provide bacterial barriers.

Biatain Ag non-adhesive foam dressings are suitable for use on fragile skin due to the absence of adhesive.

Biatain Ag adhesive foam dressings have an adhesive border based on hydrocolloid technology. Biatain Ag adhesive foam dressings shapes; sacral and heel are both specifically designed to conform to difficult to dress areas of the body. The adhesive shapes have an adhesive border based on hydrocolloid technology.

Biatain Ag foam filler/cavity is perforated to ensure flexible handling.

Biatain Ag foam dressings demonstrate in-vitro antibacterial activity for up to 7 days in bacterial strains known to be detrimental to wound healing such as *Pseudomonas aeruginosa*, *Staphylococcus aureus*, MRSA, VRE and β -haemolytic *Streptococcus*.

Sterilized by Irradiation.

Indications

Biatain Ag foam dressings are indicated for wounds with moderate to high amounts of exudate, including leg ulcers and Category II-IV pressure ulcers with delayed healing due to bacteria, or where there is a risk of infection.

Biatain Ag foam dressings may be used for second-degree burns, donor sites, postoperative wounds and skin abrasions with delayed healing due to bacteria, or where there is a risk of infection.

Biatain Ag non-adhesive foam dressings are additionally indicated for diabetic foot ulcers.

Biatain Ag foam dressings may reduce odor caused by micro-organisms in the wound. Biatain Ag foam dressings may be used to support moist wound healing on patients who are in treatment for a local or systemic infection under the discretion of a healthcare professional. Depending on the prognosis of the wound, Biatain Ag foam dressings may be used throughout the healing process to provide protection for the indicated types of wounds.

Biatain Ag foam dressings are suitable for use in combination with compression therapy.

Biatain Ag foam filler/cavity is indicated for deep exuding wounds, including leg ulcers, pressure ulcers, diabetic foot ulcers and acute wounds.

Directions

Biatain Ag foam dressings

- Rinse the wound with lukewarm water or physiological saline. Gently dry the skin around the wound. The safe use of other cleansing agents, in combination with Biatain Ag foam dressings, has not been demonstrated.
- If any film, cream, ointment or similar product is used, allow the skin to dry before dressing application.
- Remove Biatain Ag foam dressings aseptically from the packaging. Do not touch the plain (non-printed) side of the foam.
- Apply Biatain Ag foam dressings according to illustrations and application directions. Select a dressing where the foam pad overlaps the wound edge by a minimum of 2 centimeters. For small size dressings an overlap of only 1 centimeter is necessary.
- Biatain Ag foam dressings should be changed when clinically indicated or when visible signs of exudate approach the edge of the dressing. Biatain Ag non-adhesive dressings must be secured with a secondary dressing or compression bandage.
- Biatain Ag foam dressings may be left in place for up to 7 days depending on the amount

of exudate, dressing condition and type of wound.

For use of silver containing products beyond 30 days a clinical assessment by a healthcare professional is recommended.

Biatain Ag foam filler/cavity (Same directions as above point 1 and 2).

1. Apply Biatain Ag cavity according to illustrations and application directions.
2. For fixation and covering use a separate dressing, e.g. Biatain Ag adhesive foam dressings/ Biatain Ag sacral dressings (for highly exuding wounds) or Comfeel Plus Ulcer dressings/Comfeel sacral dressings or Comfeel Plus transparent dressings (for moderately exuding wounds).
3. It is recommended to change Biatain Ag cavity when the secondary dressing is changed. Ensure that Biatain Ag cavity is completely removed from the wound bed.

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Symbol

	Apply the plain (non-printed) side of the product to the wound.
WO	Work order. Production number identifying the packaging of a specific order. (LOT-number relates to production run).

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