



## Adhesive Remover Spray



Coloplast A/S  
Hoftevej 1  
3050 Humlebaek  
Denmark

### Instructions for use

23323482

Version 1 The Coloplast logo is a registered trademark of Coloplast A/S. All rights reserved.©  
2020-06-03

### Distributed by:

Coloplast Corp.  
1601 W. River Rd. N.  
Minneapolis, MN 55411 USA  
1-800-533-0464  
[www.coloplast.us](http://www.coloplast.us)

### Intended purpose

Brava Adhesive Remover spray is intended for removal of medical adhesives from intact skin.

### Cautions

Avoid direct use on non-intact skin e.g. wound beds and mucous membranes (eyes, ears, nose, mouth, vagina) as it may cause skin irritation.

Avoid direct inhalation as it may give rise to temporary irritation of the airways.

Not recommended to be used on premature infants, neonates or infants below 9 pounds / 4 kg because of their immature skin.

Extremely flammable aerosol - Danger.

Pressurized container: May burst if heated.

Keep out of reach of children.

Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

Do not spray on an open flame or other ignition source.

Do not pierce or burn, even after use.

Protect from sunlight. Do not expose to temperatures exceeding 122°F/50°C.

### Information

The product is silicone based and reusable. The product dries within seconds and does not have an effect on the adhesion of further appliances.

Ingredients: Contains a mixture of alkyl siloxanes.

Coloplast accepts no liability for any injury or loss that may arise if this product is used in a manner contrary to Coloplast's current recommendations.

### How to use

1. Spray the product around the edge of the adhesive.
2. After a few seconds gently remove the adhesive. Repeat procedure if required.

### Disposal

The product should be disposed of with other spray cans according to local regulations.

### Reporting of incidents

If, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority.

Explanation of symbols



Medical device



Indicates that the product is in compliance with European legislation for medical devices



Catalog number



Use-by date (YYYY-MM-DD)



Batch code



Date of manufacture (YYYY-MM-DD)



Manufacturer



Consult instructions for use



Indicates the upper limit of temperature to which the medical device can be safely exposed



Caution



Indicates a barcode as containing Unique Device Identification



Global Trade Item Number



Keep away from sunlight



Recyclable packaging



Compliance with Aerosol Dispenser Directive 75/324/EEC



Extremely flammable